INTENDED USE

Rapid HCV Ab Test is a chromatographic immunoassay for qualitative detection of the antibodies against hepatitis C virus (HCV Ab) in human serum or plasma samples. It is intended for use in medical institution as an aid for diagnosis and management of patients related to infection with hepatitis C as well for primary screening of blood from volunteer donors on the spot.

SUMMARY

Hepatitis C virus (HCV) is an envelope, single stranded positive sense RNA (9.5 kb) virus belonging to the family of Flaviviridae. Six major genotypes and series of subtypes of HCV have been identified. Isolated in 1989, HCV is now recognized as the major cause for transfusion associated non-A, non-B hepatitis. The disease is characterized with acute and chronic form. More than 50% of the infected individuals develop severe, life threatening chronic hepatitis with liver cirrhosis and hepatocellular carcinomas. Since the introduction in 1990 of anti-HCV screening of blood donations, the incidence of this infection in transfusion recipients has been significantly reduced. Clinical studies show that significant amount of HCV infected individuals develop antibodies to NS5 non-structural protein of the virus. For this, the third generation tests include antigens from the NS5 region of the viral genome in addition to NS3 (c200), NS4 (c200) and the Core (c22). Third generation tests have improved sensitivity and shorten the time between infection with HCV and the appearance of detectable antibodies (window period) to 60 days.

TEST PRINCIPLE

Rapid HCV Ab Test employs chromatographic lateral flow device in a cassette format. Colloidal gold conjugated goat anti-human IgM and mouse anti-human IgG are dried and immobilized on the fiberglass strip. HCV antigens are immobilized at the Test Zone (T) and goat anti mouse IgG antibodies are immobilized at the Control Zone (C). When the sample is added, it migrates by capillary diffusion rehydrating the gold conjugate. If present in sample, HCV antibodies will bind the gold conjugated anti-human IgG and/or IgM forming complexes. These complexes will continue to migrate along the strip until the Test Zone (T) zone where they are captured by the HCV antigens to form a visible red line. The un-bound gold conjugate will continue to move and bind with goat anti-mouse IgG at the Control Zone (C) forming a visible red line. If no HCV antibodies in sample, only a red line is appeared at the Control Zone (C), which indicates the validity of the test.

MATERIAL PROVIDED

1. Rapid HCV Ab Test
2. Sample buffer
3. Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED

1. Clock or timer
2. Specimen collection container
3. Centrifuge
4. Biohazard waste container

SPECIMEN COLLECTION

Fresh serum, plasma or whole blood samples can be used for this assay. Blood collected by venipuncture should be allowed to clot naturally and completely – the serum/plasma must be separated from the clot as early as possible as to avoid hemolysis of the RBC. Care should be taken to ensure that the serum samples are clear and not contaminated by microorganisms. Any visible particulate matters in the sample should be removed by centrifugation at 3000 RPM for at least 20 minutes at room temperature, or by filtration on 0.22µ filters. Plasma samples collected into EDTA, sodium citrate or heparin may be tested, but highly lipaemic, icteric, or hemolyzed samples should not be used as they could give erroneous results in the assay. Do not heat inactivate samples. This can cause deterioration of the target proteins in the sample

STORAGE

This test can be stored at room temperature (4-30°C, do not freeze!) for 18 months from the date of manufacture (see label on the pouch). The test must be used immediately after removed from the sealed pouch.


PRECAUTIONS

This test is for In Vitro Diagnostic Use only.

1. All the waste and sample should be treated in case of transmitting disease and must be properly disinfected (autoclaving is preferred) before disposal.
2. Once taking the test card out of the pouch, carry out your testing as quickly as possible (no more than 20 minutes) to avoid moisture.
3. Make sure that the test is within the validity indicated.
4. Do not modify the test procedure.
5. Avoid moisture.
6. A test giving an invalid result should be repeated.


ASSAY PROCEDURE

1. Allow the test card and sample to reach room temperature if necessary.
2. Open the pouch, Take out the test card and transfer pipet.
3. Using the transfer pipet to draw up the sample, dispense one drop (approx 40µl) specimen to the sample well marked as "S" and wait a few seconds until the sample is completely absorbed by sample pad.
4. Add one drop (approx 40 µl) sample buffer into the sample well marked as "S".
5. Read the results at 20 minutes.

Some positive samples may show positive results before 20 minutes. Results after 30 minutes may not be accurate.

INTERPRETATION OF RESULTS

POSITIVE: If two color bands are visible within 20 minutes, the test is positive and valid. The test result can be read as soon as a distinct colored band appears in test zone and control zone respectively. It indicates that HCV Ab has been detected using Rapid HCV Ab Test.

NEGATIVE: Only one red line appears in the Control Zone (C), indicating that no HCV Ab have been detected with Rapid HCV Ab Test. However, this does not exclude the possibility for infection with HCV.

INVALID: One red line should always appear at Control Zone (C) if no red line appears in the Control Zone (C), the test is invalid. The sample must be re-tested using a new device

PERFORMANCE CHARACTERISTICS

Clinical Specificity: A blood donor population of 2948 individuals was tested negative with 3 different kits from different manufacturers. The specificity of Rapid HCV Ab test was 99.55%.

Clinical Sensitivity: Among 480 clinical hepatitis C patients confirmed positive by RIBA 3.0, 479 were positive when tested with Rapid HCV Ab test. The sensitivity was 99.79%. No cross reactivity was observed with specimens from patients infected with HAV, HBV, HIV, HTLV, CMV, and TP.

LIMITATIONS

1. Negative results do not rule out the possibility of hepatitis C exposure or infection. Infection through recent exposure to HCV may not be detectable.
2. The positive result obtained with Rapid HCV Ab Test alone cannot be the final diagnosis of hepatitis C infection. Any positive result must be interpreted in conjunction with the patient clinical history and another laboratory testing results. Follow-up and supplementary testing with other analytical system (e.g. ELISA) is required to confirm any positive results.
3. This test is intended ONLY for testing of individual serum, plasma or whole blood sample. DO NOT use it for testing of other body fluids or pooled blood samples.
4. This is a qualitative assay and the results cannot be used to measure antibody concentrations.

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